ENGLISH

Intended use	For in vitro diagnostic use. For professional use only. Cystatin C Immunoparticles are intended for the quantitative determination of cystatin C in human serum, heparinized plasma and EDTA plasma by turbidimetry and nephelometry. Cystatin C measurements are used as an aid in the diagnosis and treatment of renal diseases.	
Summary and explanation	Cystatin C is a small, 13 kDa, non-glycosylated basic protein belonging to the cystatin super-family of cysteine protease inhibitors. Cystatin C is produced by virtually all nucleated cells, and is present in all investigated body fluids. The production rate is constant and is unaffected by inflammatory processes, gender, age and muscle mass (1). In the normal kidney, cystatin C is almost freely filtered through the glomerular membrane and then nearly completely reabsorbed and degraded by the proximal tubular cells. Therefore, the plasma concentration of cystatin C is almost exclusively determined by the glomerular filtration rate (GFR), making cystatin C an excellent indicator of GFR. Numerous studies and a meta-analysis incorporating 4492 subject samples have shown that serum cystatin C is superior to serum creatinine as a marker for GFR (2).	
Test principle	Human serum or plasma is mixed with the Cystatin C Immunoparticles. The resulting immune complexes are measured by turbidimetry or nephelometry. The signal generated is correlated with the concentration of cystatin C in the sample. By interpolation on a standard curve, the concentration of cystatin C in the sample is calculated.	
Reagent	Purified immunoglobulin fraction of rabbit antiserum directed against cystatin C covalently coupled to uniform polystyrene particles. Provided as a suspension and preserved with 15 mmol/L sodium azide. The Cystatin C Immunoparticle suspension is ready for use.	
Immunogen	Recombinant human cystatin C produced in <i>E. coli</i> was used as immunogen for raising the antibody coupled to the polystyrene particles.	
Precautions	<ol style="list-style-type: none">1. For in vitro diagnostic use. For professional users.2. Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.3. The possibility of antigen excess, i.e. antigen concentrations exceeding the security range of the assay, should be considered.4. This product contains sodium azide (NaN₃) a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azide. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.5. As with any product derived from biological sources, proper handling procedures should be used.	
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the user must verify the conditions. If unexpected results are obtained which cannot be explained by variations in laboratory procedures and a problem with the reagent is suspected, contact our Technical Services. For on-board stability data please refer to the Specific Performance Characteristics for the specific instrument.	
Specimens	<p>Sample material: Human serum, heparinized plasma and EDTA plasma.</p> <p>Samples can be stored at 2-8 °C, if analysed within 2 days. Otherwise, samples must be stored at -20 °C. Samples are stable for at least 6 months at -20 °C. It is recommended that frozen samples are thawed at 37 °C; samples must be mixed well before analysis.</p> <p>As with any product derived from biological sources, proper handling procedures should be used.</p> <p>For interference please refer to the Specific Performance Characteristics for the specific instrument.</p>	
Procedure	<div><div><u>Materials provided</u> Dako Cystatin C Immunoparticles</div><div>Code No. LX002</div></div> <div><div><u>Materials required but not provided</u> Dako Cystatin C Calibrator Dako Cystatin C Control Set Dako Reaction Buffer 9 Diluent</div><div>Code No. X0974 Code No. X0973 Code No. S2361 See the Application Note for the individual instrument</div></div> <div><u>Assay procedure</u> Detailed, validated Application Notes for the determination of cystatin C in human sample material, based on optimized test systems (3), are available from Dako for a range of automated instruments. Applications Notes for Hitachi 917 and IMMAGE® Immunochemistry System are included on the following pages. For instruments not covered by a specific Application Note, a General Application Note, which provides general assay procedure guidelines is available. When using the General Application Note, each individual laboratory should validate their test system. Please consult the instruction manual of the instrument used regarding precautions, operation, maintenance, and instrument performance characteristics.</div> <div><u>Establishment of the calibration curve</u> The assigned value of the Dako Cystatin C Calibrator, code No. X0974, is given on the Analytical Value Sheet included with the calibrator. The exact measuring range depends on the specific lot of the calibrator and the Application Note for the individual instrument.</div>	

Calibration curves are constructed by multi-point calibration. Depending on the instrument, a dilution series is either prepared automatically by the instrument or manually, and then applied on the instrument. For details on interval for recalibration, please refer to the following pages.

Assay of samples

Samples are measured according to the Application Note for the individual instrument.

If the results obtained are outside the measuring range, the assay can be repeated using a post concentration or post dilution of the sample as stated in the Application Note.

Internal quality control

Dako Cystatin C Control 1 and Control 2, code No. X0973, must be assayed with each establishment of a calibration curve and with each run of samples. The controls should be assayed and evaluated as samples. The assigned values and confidence intervals are given on the Analytical Value Sheet included with the controls. If quality control results fall outside the range stated on the Analytical Value Sheet or outside the range established by the individual laboratory, recalibration is recommended.

Results The analytical results are automatically calculated by the instrument.

Limitations For interference please refer to the Specific Performance Characteristics for the specific instrument.

Expected values Individuals 1-50 years: 0.55-1.15 mg/L
Individuals >50 years: 0.63-1.44 mg/L

The reference intervals were determined from a population of 69 subjects ≤ 50 years and 94 subjects > 50 years of age, all with normal glomerular filtration rate. The reference intervals were calculated using non-parametrical statistics and represent the central 95% of the population. Healthy individuals ≤ 50 years of age have generally a glomerular filtration rate above 80 mL/min/1.73 m², whereas the glomerular filtration rate of healthy individuals > 50 years decreases with age (4). The above reference intervals may be used as a guide for demonstrating normal or reduced GFR, but it is recommended to determine the reference intervals for the local population.

Specific performance The Dako Cystatin C Assay is based on the usage of Dako Cystatin C Immunoparticles **characteristics** together with Dako Cystatin C Calibrator, Dako Cystatin C Control Set, Dako Reaction Buffer 9, and the diluent recommended in the Application Note for the specific instrument. Please refer to the following pages for the Specific Performance Characteristics for the Dako Cystatin C Assay on Hitachi 917 and IMMAGE® Immunochemistry System.

Specificity The specificity of the cystatin C antibody, before coupling to the particles, was ascertained by crossed immunoelectrophoresis. Only the cystatin C arch appeared when the antibody was tested against concentrated human urine from patients with tubular proteinuria.

Method comparison Results obtained with the Dako Cystatin C assay were compared with a commercially available nephelometric assay.


Estimation of GFR from cystatin C Measurement of the serum/plasma cystatin C concentration allows for estimation of GFR (5). Plasma cystatin C concentration has been measured with the Dako Cystatin C assay on 536 patients, who underwent iothexol clearance determination (4). The patients covered an age range from four months old to 93 years old.

The following equation for estimation of GFR has been established (5):

$$\text{GFR (mL/min/1.73 m}^2\text{)} = [84.69 \times \text{cystatin C (mg/L)}^{-1.680}]^*$$

The equation results in the following relationship between cystatin C and GFR:

Cystatin C (mg/L)	GFR (mL/min/1.73 m ²)*
0.6	200
0.7	154
0.8	123
0.9	101
1.0	85
1.1	72
1.2	62
1.3	55
1.4	48
1.5 – 1.6	41
1.7 – 1.8	33
1.9 – 2.0	28
2.1 – 2.3	23
2.4 – 2.6	18
2.7 – 3.0	15
3.1 – 4.0	10

* If the patient is < 14 years of age, the result has to be multiplied with a prepubertal factor of 1.384.

The equation and table above provide estimates of GFR adjusted to a standard body surface area of 1.73 m² (mL/min/1.73 m²). To obtain estimates of absolute GFR (mL/min), multiply the result obtained in the equation or table above with the patient's estimated body surface area (9) divided by 1.73. The equation and table may be used as guidelines for estimating GFR, but it is recommended that each individual laboratory verify the equation.

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
1. Grubb AO. Cystatin C – properties and use as diagnostic marker. Adv Clin Chem 2000;35:63-99.

2. Dharnidharka VR, Kwon C, Stevens G. Serum cystatin C is superior to serum creatinine as a marker of kidney function: A meta-analysis. Am J Kidney Dis 2002;40:221-6.
3. Blirup-Jensen S. Protein standardization III: Method optimization. Basic principles for quantitative determination of human serum proteins on automated instruments based on turbidimetry or nephelometry. Clin Chem Lab Med 2001;39:1098-1109.
4. Nilsson-Ehle P. Iohexol clearance for the determination of glomerular filtration rate: 15 years' experience in clinical practice. eJIFCC vol 13 no2: <http://www.ifcc.org/ejifcc.htm>.
5. Grubb A, Nyman U, Björk J, Lindström V, Rippe B, Sterner G, et al. Simple cystatin C-based prediction equations for glomerular filtration rate compared with the modification of diet in renal disease prediction equation for adults and the Schwartz and the Counahan-Barratt prediction equations for children. Clinical Chemistry 2005;51:1420-31.
6. Evaluation of precision performance of clinical chemistry devices. Approved guideline. NCCLS Document EP5-A, 1998.
7. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-85.
8. EURACHEM/CITAC working group. EURACHEM/CITAC Guide: Quantifying uncertainty in analytical measurement. Eurachem CITAC, second edition 2000; P1:08-111.
9. DuBois D, DuBois EF. A formula to estimate the approximate surface area if height and weight be known. Arch Int Med 1916;17:863-71.

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Hitachi 917

Specific Performance Characteristics for Dako Cystatin C Assay

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Application Note for Determination of Human Cystatin C
Instrument Settings for Hitachi 917 (Software Version 7146020-02-48)

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Diluent – Saline

Reaction buffer (R1) – Dako Reaction Buffer 9, code No. S2361

Immunoparticles (R3) – Dako Cystatin C Immunoparticles, code No. LX002

[*] Defined by the customer.

[**] The calculated relative concentration of Std.1 and Std.6. In order to get a warning for samples with concentration below the concentration of the lowest standard, the lower technical limit has to be 0.001 mg/L higher than the lowest standard.

- [***] The concentration is calculated as the factor times the calibrator value for the specific lot of Dako Cystatin C Calibrator, code No. X0974, (C(cal), stated in mg/L on the X0974 Analytical Value Sheet). Number of decimals stated for Std.1 defines the number of decimals in printout. Two decimals are recommended.

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IMMAGE® is a registered trademark of Beckman Coulter, Inc. Fullerton, CA. 92835.

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Application Note for Determination of Human Cystatin C

Instrument Settings for IMMAGE® Immunochemistry System and the IMMAGE® 800 Immunochemistry System

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Procedure

1. Program a user-defined reagent with the parameters listed below. For more detailed instructions, please refer to chapter 5 of the IMMAGE® or IMMAGE® 800 Immunochemistry System Operations Manual.
2. Transfer Dako Cystatin C Immunoparticles to compartment "A" of a new User-Defined Cartridge.
3. Transfer Dako Reaction Buffer 9 to compartment "B" of the cartridge. Please see the note below.
4. Name the Chemistry "CYST".
5. The calibrator value is lot dependent. When using the Dako Cystatin C Calibrator (code No. X0974) the following dilution series is prepared manually according to the dilution scheme shown below.
6. Use IMMAGE® DIL1 (Beckman Coulter PN 447255) as the sample diluent and configure it on the IMMAGE® as DIL10.
7. For the IMMAGE® 800 the sample diluent (Beckman Coulter PN 447255) may be configured as DIL1.

Dilution scheme

Standard No. (i)	X0974 μL	NaCl 0.9% μL	$f_{(i)}$	Conc. Factor	$C_{(i)}$ Concentration
Std. (1)	[20]	[360]	$f_{(1)}$	0.0526	[$0.0526 \times C_{(cal)}$]
Std. (2)	[20]	[180]	$f_{(2)}$	0.1000	[$0.1000 \times C_{(cal)}$]
Std. (3)	[40]	[120]	$f_{(3)}$	0.2500	[$0.2500 \times C_{(cal)}$]
Std. (4)	[75]	[100]	$f_{(4)}$	0.4286	[$0.4286 \times C_{(cal)}$]
Std. (5)	[90]	[60]	$f_{(5)}$	0.6000	[$0.6000 \times C_{(cal)}$]
Std. (6)	[100]	[0]	$f_{(6)}$	1.0000	[$1.0000 \times C_{(cal)}$]

The concentration of cystatin C in the standards ($C_{(i)}$) (mg/L) is calculated from the cystatin C concentration ($C_{(cal)}$) (mg/L) stated for the specific lot of the calibrator (X0974) and the relative concentration factor ($f_{(i)}$) by multiplication, i.e. $C_{(i)} = f_{(i)} \times C_{(cal)}$. Please see the Analytical Value Sheet for the exact $C_{(cal)}$ value of the current lot of Dako Calibrator, code No X0974.

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Chem. Name		Units	mg/L
Lot Number	XXXXXXXX	Protocol	Non-Competitive Nephelometric
Reagent Serial Number	XXXX	Reagent Expiration Date	xx/xx/xx
Sample or Dilution Volume	10 μL	Gain	1
Reaction Buffer Volume	0 μL	Cal. Dilution	1:5
Compartment A Volume	33 μL	Sample Dilution	1:5
Compartment B Volume	195 μL	Reaction Time	1.5 minutes
Calibrator Levels/Replicates	6/2	Calibrator Set Points	Lot-specific entries

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Note

The maximum volume of compartment "A" and "B" of the UDR cartridge is 8.0 mL, and it has a dead volume of 250 μL. While the system dispenses 195 μL of buffer from compartment B according to the protocol described above, it actually aspirates a total volume of 205 μL to avoid any possibilities of liquid contamination. Thus when the cartridge is maximally filled, approximately 35 tests can be run. It is advisable to monitor the reagent usage, and refill compartments A and B when appropriate until the full 300 tests from the cartridge have been utilized.

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IMMAGE[®] 800 protocol information

Information Protocole IMMAGE[®] 800

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Chem. Name		Units	mg/L
Reagent Lot Number	XXXXXXXX	Protocol	Non-Competitive Nephelometric
Cartridge Lot Number	XXXXXXXX	Cartridge Serial Number	XXXX
Tests per Cartridge	35 tests	Reagent Expiration Date	XX/XX/XX
Sample or Dilution Volume	10 µL	Gain	1
Reaction Buffer Volume	0 µL	Cal. Dilution	1:5
Compartment A Volume	33 µL	Sample Dilution	1:5
Compartment B Volume	195 µL	Reaction Time	1.5 minutes
Update Level	4	Replicates	2
Calibrator Levels/Replicates	6/2	Calibrator Set Points	Lot-specific entries

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